



YZ20P5
OPERATION MICROSCOPE
INSTRUCTION MANUAL

66 VISION TECH CO., LTD.

THE PEOPLE'S REPUBLIC OF CHINA

Preface

Thank you for purchasing our YZ20P5 Operation Microscope. Please read this specification carefully for the sake of your better use.

General Requirement for Safety

Please read carefully about following precautions to avoid unexpected personal injury as well as the product being damaged and other possible dangers.

Precautions

1. Do not use this instrument in the environment where is prone to fire and to blast or where there is much dust and high temperature. Use it in the room and simultaneously be careful to keep it clean and dry.
2. Check that all the wires are correctly and firmly connected before use. Ensure that the instrument is well grounded.
3. Please pay attention to all the rating of the electrical connecting terminal.
4. Please only use fuse according to the specifications and rated stipulated by our product.
5. Use the power cable supplied with this instrument.
6. Don't touch the surface of the lens and prism with hand or hard objects.
7. Turn off the main power first before replacing the main bulb and fuse.
8. To prevent the instrument from falling down to floor, it should be placed on the floor that is horizontal.
9. Turn off the power and cover the instrument with dustproof cover when it is not in use.
10. In case there is any trouble, please first refer to the trouble-shooting guide. If it still can't work, please make contact with the authorized distributor or our Repair Department.

THE SAFETY MARKS USED IN THIS INSTRUMENT



TYPE B



ATTENTION



TERMINAL OF
THE PROTECTIVE
GROUNDING



SEPARATE
COLLECTION



UPWARD



CAUTION
PLACE



AVIOD DAMP



SCALD WARNING

Contents

1. Features & Specifications	1
1.1 Features	1
1.2 Specifications	1
2. The name and use of components	3
3. Assembly	5
3.1 Assembly of stand pillar	6
3.2 Assembly of light source box	7
3.3 Assembly of microscope	7
3.4 Assembly of manipulating handle	7
3.5 Assembly of fiber optics	7
4. Preparations for use	7
5. Use of instrument	8
5.1 Use of foot control switch (Fig. 3) (You can choose manual focus.)	8
5.2 Installation & adjustment before the use	8
5.3 Adjustment during using	9
5.4 Movement & storage after the using	10
6. Maintenance	10
6.1 Replacing the consumable parts	10
6.1.1 Replacing the bulbs	10
6.1.2 Replacing the fuse	11
6.2 Cleaning & sterilization	11
6.3 Trouble-shooting guide	12
6.4 Order consumable parts	12
7. Responsibility	13
8. Transportation, Storage & Rejection	13
9. Spare parts & tools of the instrument	13
10. Assembly of the optional accessories (additional order)	14
10.1 Erect gminate ocular drawtubes (for ENT)	14
10.2 Beam splitter, CCD adapter	14
10.3 10x ocular, 16x ocular	15
10.4 Universal coupling (for ENT)	15
10.5 Objective, F100mm/F175mm/F250mm/F300mm/F350mm/F400mm	16
10.6 Clamp rack, fission column	17
10.7 Handy casing	17
10.8 lengthen arm	17

1. Features & Specifications

1.1 Features

The YZ20P5 Operation Microscope is a portable microscope. Its small and portable machine body makes it highly agile. It is sure to meet the requirements for general ophthalmic micro-surgical.

The microscope is a single binocular one, with three magnification steps (standard assemblage), offering sharp image and large field of vision. A cold light source is used in the illumination system, which is harmless to the tissue of eye. With quality halogen bulbs, PHILIPS or OSRAM, being used, the illumination is bright and even. An intensifying red-reflex module is added.

A spring balance system is designed for the arm so that the microscope can move upwards and downwards stopping at any desired position. Fine focusing adjustment can be performed. The adjusting functions of the equipment include magnification, focusing and inclination, while focusing can also be controlled by foot switch.

The environment requirements for using this equipment is temperature +10°C to +40°C, relative humidity 30% to 75%, atmospheric pressure 50hPa to 1060hPa.

1.2 Specifications

(1) Parameter of microscope

Focus length of objective $f = 200\text{mm}$

Eyepieces $12.5\times/16\text{B}$

Magnification & Visual Field:

Objective focus	Five magnification steps parameter					Illumination field (mm)
	Three magnification steps parameter					
magnification (Field diameter)mm	3.2X (ϕ 63)	5.3X (ϕ 38)	8X (ϕ 25)	12X (ϕ 17)	20X (ϕ 10)	ϕ 37.5mm

(2) Working distance 190mm

(3) Parameter of Binocular tube

Inclined binocular tubes 45°

Diopter adjustment range $\pm 5\text{D}$

Pupil distance range $50\text{mm}\sim 70\text{mm}$

Height of eyecups 18mm

(4) Parameter of Illumination

6°+0° coaxial illumination

Maximum intensity 30000Lx or up

(5) Parameter of Position adjustment

Maximum stretching radius of arm 870mm

Vertical movement range (from floor to front surface of big objective) 700mm~1100mm

Range of fine focusing adjustment 30mm

(6) Parameter of Electric

Input voltage AC220V ± 10%/50Hz, or
AC110V ± 10%/60Hz

Input power 120VA

Fuse T type AC250V 1.25A、 T type AC125V 2.5A
(For example 51S-013H、 51S-025L)

Bulb 12V/100W
(Medical halogen lamp Bulb with cold reflection)

Electrical safety Conform to standard IEC601-1 Type I, Class B

2. The name and use of components




Fig. 1

[1] Caster

Move and support the equipment

[2] Foot control switch

Be used to control the focusing functions of microscope

- [3] Floor stand
Be used to support & fix stand pillar
- [4] Stand pillar
- [5] Fuse
Two (1.25A/2.5A) are used to power supply
- [6] Light-adjusting knob
Used for continuous adjustment of illumination intensity.
- [7] Power supply switch
Be used to switch on or switch off power supply.
- [8] Spare bulb module
In case the working bulb burnt during an operation, the spare lamp module may be inserted so as to ensure the operation to be continued smoothly.
- [9] Bulb module
- [10] 7-pin plug/socket
Connecting wires for focusing control
- [11] Fixing nut
This nut lets microscope suspend on the small horizontal arm. Though there is protective melt, you must check this nut tightened or not.
- [12] Star fixation knob (with sterilized cap)
Fix the revolving angle of the suspended spindle of microscope.
- [13] Star fixation knob (with sterilized cap)
Lock the microscope and make it not revolve at up and down plane.
- [14] Diopter adjusting ring
Adjust ocular diopter by rotating this ring. The range of adjustment is $\pm 6D$.
- [15] Eyecups
Adjust exit-pupil distance. Its height is 18mm and it may be taken off or rolled down.
- [16] Manipulating handle (with sterilized cap)
For rough focusing, move the microscope up and down or right and left.
- [17] Magnification knob (with sterilized cap)
With different objectives the magnification factor of each step differs. Rotate the knob to change the magnification.

[18] Fixing block circle

Fix the manipulating handle so that it will not loose.

[19] Fuse pin

Avoid the microscope dropping when the fixing nut looses while unloading it or suspending it.

[20] Fiber optic

Lead the light beam from bulbs to the operating position.

[21] Star fixation knob(with sterilized cap)

Lock the little arm by tightening it firmly so as to prevent the microscope from moving up and down.

[22] Star fixation knob

Tighten it in order to lock the little arm and make the microscope not move vertically.

[23] 5-pin plug/socket

Connect the foot switch.

[24] Star fixation knob

Tighten it in order to lock the light source box to make it not revolve when moving or storing.

[25] Inner hexagonal tightening bolt

Lock the stand pillar and the light source box

[26] 110/220V selecting switch.

[27] Power plug

Input power supply

3. Assembly

This equipment may be installed either by user self with reference to the manual or by the servicemen sent by the manufacturer or authorized representatives when facing real difficulties.

This equipment is packed in one package. Please open the packages in the direction indicated by the mark on the packages. Take out all the parts and assemble them according to the following procedures:

3.1 Assembly of stand pillar




Fig. 2

Please take out the floor stand [3] from the package, lay it on the ground.

Take out the stand pillar [4], turn out the inner hexagonal bolt and gasket of its end, insert it into the hole of the floor stand [3], and then turn the stand pillar [4] to make the column pin on the base support clip in the groove on the end of stand pillar, assemble in order the unassembled gasket、spring gasket、inner hexagonal blot, and fasten it firmly with 10mm hexagonal wrench (Fig.2).

3.2 Assembly of light source box

Take out the light source box while not loosen the tightening bolt of the star knob [21] in order to avoid the small arm rebounding and harm people, set the bottom of the light source box into the stand pillar [4] and insert the inner hexagonal tightening blot [25] into the groove of the end of stand pillar [4] and fasten it firmly with 6mm hexagonal wrench(Fig.2).

3.3 Assembly of microscope

- (1) First check the tightening blot of star knob [12], make sure that the end of the bolt does not protrude out of the hole surface.
- (2) Remove the fixing nut [11]. Pull out the safety pin [19] with one hand, hold the connecting part of microscope with another hand. Insert the shaft into the hole of the arm from the bottom to top. Then release the safety pin [19], it will spring back to the groove on the shaft. Finally fasten the fixing nut [18] (Fig. 2).
- (3) Insert the 7-pin plug [10] into the 7-pin socket. Pay attention to align marks for alignment on both the plug and socket.

3.4 Assembly of manipulating handle

Insert the manipulating handle [16] into the holes on the two sides of the microscope holder respectively. The proper handle position should be down to an angle of 40° (as shown in Fig. 1). Then fasten the fixing block circle on the project plane of the two sides of the microscope and make the manipulating handle not to move axile direction (Fig. 2).

3.5 Assembly of fiber optics

- (1) Remove the fiber optics protective cover.
- (2) Insert one end of the fiber optics [20] into the side hole of the light source box and insert another end into the hole behind the microscope (Fig. 2).

4. Preparations for use

- (1) Check whether the mains voltage, frequency complies with what required by the equipment. If not, do not start it.



Caution: Set the input voltage at 220V. When it is 110V, please switch the selecting switch [26] to the 110 block and change the fuse[5] as type AC125V 2.5A fuse tube we supplied.

- (2) Check the grounding of power supply. Make sure the equipment have a good ground-wire connection.
- (3) The equipment comes with a three-core power cable. Please select matched

power socket.



Caution: Please use the power cable provided by the manufacturer or the power cable according with IEC227 standard to ensure that the equipment is well grounded.

- (4) When the powder switch [7] on the power controlling box is pressed at "ON" position, power is on. When pressed at "OFF" position, power is off. The switch should always be "OFF" before the power wire is connected with power socket.
- (5) Insert the plug of the power wire of the equipment into the mains' outlet (It must be well grounded).
- (6) Switch on the power. Check whether the illumination system works properly.
- (7) After checking, please turn off power switch on power controlling box and disconnect the power plug and cover the equipment with dust cover.

5. Use of instrument

5.1 Use of foot control switch (Fig. 3) (You can choose manual focus.)




Fig. 3

- (1) First connect the 5-pin plug on the foot control switch [2] with the 5-pin socket [23] under the power controlling box (Fig.1). Pay attention to align marks for alignment on both the plug and socket. Then connect power cable with mains outlet and switch on the power switch [7]. so that you start foot control.
- (2) For fine focusing, step down the left pedal of the foot control switch, the microscope moves down; while step down the right pedal, the microscope moves up (Fig.3).

5.2 Installation & adjustment before the use

- (1) Sterilize the manipulating handles and sterilizable caps in advance.
- (2) Make the little arm to be in a horizontal position. Adjusting the height, ensure that the objective is about 200mm from the surgical site.
- (3) Switch on power supply, check the bulbs. If any bulb burns, replace it

immediately.



Caution: Burnt bulbs must be replaced before an operation to ensure that the operation may be performed without interruption.

5.3 Adjustment during using

- (1) First move the equipment to a proper position and skid the brake to prevent the equipment from moving easily. Connect the power plug with power socket [27] on the microscope.
- (2) Place the foot switch [2] to a suitable position. Connect the 5-pin plug of the foot switch with the socket under the power controlling box.
- (3) Insert the plug of power cable under the power controlling box into the mains outlet. Then switch on power switch [7].
- (4) Move the fine focusing drag to the starting position, where the fine focusing indicator point to the central point of the drive, by using the foot switch [2] or the focusing knob.
- (5) Illumination intensity adjustment. Rotate the light adjusting knob [6] on the power controlling box clockwise to increase the intensity, counterclockwise to reduce the intensity. Adjust it until appeal the operation.
- (6) Rough focusing. Loose the star fixation knob [21], hold the handle [16] to move the microscope upwards or downwards to position the light spot on the surgical site and adjust the focus till the image is clear with 8x magnification. Manipulating method referes to 5.1(2).
- (7) Eyecup adjustment. The eyecup [15] can reduce the observing disturbance coming from the foreign light. The height is limited to 18mm. If a doctor wearing glasses makes the operation, need to roll down the eyecup, that is turn the rubber eyecup over the ocular tube.
- (8) Diopter adjustment. The adjustment range of the diopter adjusting ring [14] of the ocular is $\pm 6D$, with 1D for each division. Rotate the diopter adjusting ring, match the value on the scale of the diopter adjusting ring with the diopter of the surgeon. If the surgeon wears glasses, align "0" position on the diopter adjusting ring to the white line mark on the ocular tube because glasses has corrected the surgeon's vision.
- (9) Pupil distance adjustment. While adjusting the pupillary distance, you can observe while rotate the pupillary distance adjusting knob, until both eyes see properly, visual fields coincide and stereoscopic effect obtained. If the pupillary distance of the surgeon is known, adjust directly to this value indicated on the surface of the PD adjusting rod.

★The above steps should be done before surgeon's personal sterilization.

- (10) Move the microscope upwards. Put the sterilizable caps at the following places:

Star fixation knob [12], [13], [21] Manipulating handle[16] Magnification knob [17]

★ Disposable sterile cover may be used to cover the microscope at the surgeon's choice. After that, move the microscope to position the light spot on the center of the surgical site.

5.4 Movement & storage after the using

- (1) Put off all of the sterilized cap and knob and sterilize them for the next use.
- (2) Draw the microscope back to the nearest position to the stand. Fasten every star knob firmly so as to fix the arm and the microscope.
- (3) Disconnect power wire and foot switch plug, store in a suitable place.
- (4) Loose the brake of the caster[1] before moving.
- (5) Grasp the power control box firmly with two hands in order to move the equipment slowly and carefully meanwhile avoid bumping and leaning.
- (6) Replace burnt bulb if any for the next use. Refer the method to 6.1.1.

6. Maintenance

6.1 Replacing the consumable parts.



Caution: Replaced parts are treated as ordinary industrial rubbish.

6.1.1 Replacing the bulbs

- (1) The lamp house is equipped with a spare bulb module [8]. In case the bulb burnt out during an operation, only need to turn off power switch and pull out the bulb module [9](Fig. 7) then insert the spare bulb module [8] to the same direction as being pulled out. After that, switch power on and continue the operation.




Fig. 4

- (2) After the operation, replace the burnt bulb immediately for next operation. While replacing, pull out ceramic socket and take the bulb out of two clamping springs. Then replace a new bulb in a reverse way (Fig. 4).



Caution: watch out for hot bulbs.

- (3) Insert the bulb module to the same direction it was removed.



Caution: Please use the special bulb for the equipment

6.1.2 Replacing the fuse

Remove the fuse holder under the power controlling box. Replace it with a new fuse, then fasten the cover.



Caution: Only use fuses of the same type, specification and rated value. To be safe, please disconnect power from the mains outlet before remove the lower cover and replace the fuse.

6.2 Cleaning & sterilization

- (1) Every equipment has been fully checked before delivery ensuring its proper performance. However, a proper maintenance is necessary. It should not be disassembled by unskilled or unqualified technicians. Otherwise the equipment may be damaged and the performance may be affected. It is advised by the manufacturer that the equipment takes regular maintenance at three month's interval or according to the specific condition.
- (2) Do not place the equipment in a dusty, moist or corrosive environment to avoid damage to the equipment.
- (3) Do not disassemble lenses yourself. If there is dust stained on the lens, blow them with a blowball or brush them with a dust pen. Greasy or water stains can be cleaned off with mirror-cleaning paper or a drop of liquid solvent (1:1 mixture of C_2H_5OH and CH_3OCH_3), then blow it carefully. Be careful to prevent the solvent from infiltrating the edges of lens.

- (4) The outer surface of the equipment may be cleaned with wet cloth. The remaining stains can be cleaned off with the mixture of 50% C₂H₅OH and 50% H₂O (distilled). Please not use corrosive cleaning agent.
- (5) Eyecups should be cleaned with water. Place them on the eyepieces after dry. All sterilizable caps should be autoclave sterilized. Suggested conditions are as following: temperature: 120℃, time: 20min; or temperature:134℃, time: 5min.
- (6) The accessories not in use should be put in a closed box with desiccant.

6.3 Trouble-shooting guide

In case there is any trouble, please check according to the following table. If it still does not work, please contact the Repair Department of 66 Vision Tech Co., Ltd.

Table 1 Trouble Shooting Guide

Trouble	Possible Cause	Remedy
Illumination doesn't light	The cable isn't connected correctly with the power socket.	Connect the power cable correctly.
	The position of the light adjusting knob is wrong after switching on power switches.	Adjust the light-adjusting knob.
	The bulb plug doesn't connect correctly with its socket.	Take out bulb, scrape oxide. Then reinstall, screw the locking bolt firmly.
	The bulb has burnt (The mains exceed voltage rating).	Replace new bulb (stabilize the voltage rating).
	The fuse has blown.	Change the fuse.
The foot switch does not work.	The plug of foot switch is not well connected with the socket under the powder controlling box. The two connectors on the second arm are not well connected.	Connect again.
The light spot is too dark or not even.	The bulb socket is bad plugged.	Plug it correctly.
	The fiber optics is not inserted to the end.	Insert it to the end.

6.4 Order consumable parts

Name	Specification	Suggested Type
Medical halogen bulb with cool reflection.	AC12V/100W	HLX64627 (OSRAM)
Fuse	T type AC250V 1.25A T type AC125V 2.5A	51S-025L,51S-013H (HOLLYLAND CO.,LTD)

7. Responsibility

We may provide the circuit diagram, electric component list and other details of the equipment at the request of customers to meet their need for repairing.

If any information, service or consultation is needed, please contact with us directly or with our authorized distributors.

8. Transportation, Storage & Rejection

During the transportation, take care to protect it from wetness, upside down and violent vibration. The relative humidity should be 10% to 90%, environmental temperature -10°C to +40°C and atmospheric pressure 500hPa to 1060hPa.

This equipment should be stored in a well-ventilated room without corrosive gas where the relative humidity should be 10% to 90% and environmental temperature -10°C to +40°C. If the assembled equipment should be moved or transported in short distance, please lock all the movable parts (refer to 5.4). The inclination angle should be less than 10°. If for long distance transportation, please repack it with original packages.

If it is stored for over five years, please contact with us directly or with our authorized distributor to check the equipment again.

Discard the equipment according to the local environment protection regulations. Do not pollute environment.

9. Spare parts & tools of the instrument

1. Sterilized cap for Star fixation knob	4
2. Sterilized cap for magnification knob	2
3. Sterilized cap for Manipulating handle	2
4. Medical halogen bulb with cool reflection (12v 100w)	1
5. T type AC250V 1.25A; T type Ac 125V 2.5A fuse tube	4/each
6. 6mm inner hexagonal wrench	1
7. 10mm inner hexagonal wrench	1
8. P54M bulb socket	2

10. Assembly of the optional accessories (additional order)

10.1 Erect geminate ocular drawtubes (for ENT)

Unscrew the inner hexagonal tightening bolt [A3] by the 2mm inner hexagram screw drive. Take down the 45° oblique geminate ocular drawtubes [A2] and fix the erect geminate ocular drawtubes [A1]. Then screw down the inner hexagonal tightening bolt [A3] by the 2mm inner hexagram screw drive. (Fig.5)




Fig. 5

- [A1] Erect geminate ocular drawtubes
- [A2] 45° oblique geminate ocular drawtubes
- [A3] Inner hexagonal tightening bolt
- [A4] Objective




Fig. 6

- [B1] beam splitter
- [B2] screw flange
- [B3] CCD adapter
- [B4] focus adjusting ring
- [B5] screw flange
- [B6] CCD connector
- [B7] dust-proof cover
- [B8] CCD
- [B9] dust-proof cover
- [B10] dust-proof cover

10.2 Beam splitter, CCD adapter

- (1) Unscrew the inner hexagonal tightening bolt [A3] by the 2mm inner hexagram screw drive.
- (2) Fix the beam splitter [B1] in the microscope.
- (3) Select one beam split-port in the left or right side of the beam splitter as the position of the CCD adapter [B3]. Unscrew the dust-proof cover [B10] on the beam split-port and keep it well.
- (4) Take out the CCD adapter [B3] and screw off the dust-proof cover [B7] [B9]. Fix the CCD adapter [B3] into the microscope beam split-port and wedge the orientation bolt of the beam split-port into the orientation groove of the CCD adapter [B3]. (There are three grooves on the plane. Select one according to your need.) Screw down the screw flange [B2] on the top of the adapter.
- (5) Screw off the lens and screw flange of the CCD [B8] and keep them well.

Screw the CCD connector [B6] into the CCD lens and tighten it.

- (6) Screw down the CCD with the CCD adapter [B3].
- (7) Connect the power plug to the appointed power.
- (8) Connect one end of the video frequency wire to the CCD signal fan-out. The other end should be connected to the Video In port of the monitor or video camera. Switch on the power after all have been done. Observe the image on the monitor and rotate the focus-adjusting ring [B4] to make the image clear. Attention to the orientation of the image. Rotate the CCD out by inches till the image's orientation is right if it is declining. Screw down the screw flange [B5] to lock it.
- (9) Fix the geminate ocular drawtubes to the beam splitter [B1] and screw the inner hexagonal tightening bolt [A3] by the 2mm inner hexagram screw drive. (Fig.6)



Attention: After the above assembly, the dynamic balance of the small arm should be adjusted. Deploy the small arm horizontally, and drag the microscope up and down, and then compare the resistance between upwards moving and downwards moving. Rotate the balance adjustment bolt counter clockwise by the single-groove screw tool [C1]. If the resistance of downwards moving is bigger, then rotate the balance adjustment bolt clockwise by the single-groove screw tool. Do this repeatedly till the resistances of upwards moving and downwards moving are almost equivalent. Then the small arm may stay steadily in any position of $\pm 20^\circ$ relative to horizontal line. (Fig.8)

10.3 10x ocular,16x ocular

Pull out the 12.5x ocular from the geminate ocular drawtubes and insert the 10x ocular, 16x ocular. The magnification is as 0.8 times,1.25 times as before. The size of visual field is as the times as before.



Attention: After the above assembly, the dynamic balance of the small arm should be adjusted. Deploy the small arm horizontally, and drag the microscope up and down, and then compare the resistance moving are almost equivalent. Then the small arm may stay steadily in

10.4 Universal coupling (for ENT)

It is used for ENT that the universal coupling [C1] and the erect geminate ocular drawtubes' [A1] combination. The microscope lens can be rotated in X, Y, Z plane after the star fixation knob [C2] is unscrewed. (Fig.7)




Fig. 7

- [C1] Universal coupling
- [C2] Star fixation knob




Fig. 8

- [D1] Single-groove screw tool
(provide for yourself)

10.5 Objective,

F100mm/F175mm/F250mm/F300mm/F350mm/F400mm

Unscrew the objective [A4] under microscope and fix the special objective [A4].
(Fig.5)

The size of visual field according to next table:

Objective	Five magnification steps parameter					Illumination field (Work Distance)
	Three magnification steps parameter					
$f'=100\text{mm}$	6.4X (φ 31.5)	10.6X (φ 19)	16X (φ 12.5)	24X (φ 8.5)	40X (φ 5)	φ 19mm (90mm)
$f'=175\text{mm}$	3.6X (φ 55)	6X (φ 33)	9X (φ 22)	13X (φ 15)	23X (φ 9)	φ 33mm (170mm)
$f'=250\text{mm}$	2.6X (φ 79)	4.2X (φ 47.5)	6.4X (φ 31)	9.6X (φ 21)	16X (φ 12.5)	φ 47mm (240mm)
$f'=300\text{mm}$	2X (φ 94.5)	3.5X (φ 57)	5.4X (φ 38)	8X (φ 25)	13.3X (φ 15)	φ 57mm (290mm)
$f'=350\text{mm}$	1.8X (φ 110)	3X (φ 67)	4.6X (φ 44)	6.8X (φ 30)	11.5X (φ 17.5)	φ 67mm (340mm)
$f'=400\text{mm}$	1.6X (φ 126)	2.6X (φ 76)	4X (φ 50)	6X (φ 34)	10X (φ 20)	φ 74mm (390mm)

10.6 Clamp rack, fission column

Fix the clamp rack [E2] on one side of the table [E3]. Make the inboard side of the clamp rack [E2] cling to the edge of the table before clamp. Screw the lead screw to make the cushion to be pressed under the table. Lift the spanner 90° relative to the lead screw and clamp them. Assembly of the lamp-house box [E1] is indicated in 3.2. (Fig.9)

Assembly of the underside of the column [F3] is indicated in 3.1. Fix it to the base [3]. Screw the upside of the column [F2] into the underside of the column. Screw down the upside through the spanner hole [F1]. (Fig.11)

10.7 Handy casing

The handy casing is composed of three canvas-bags, a luggage barrow (ordering them according to your need).

10.8 lengthen arm

Lengthen arm can be used to increase stretch radius.




Fig.9

[G1] Star fixation knob (with sterilized cap)

Lock the second arm by tightening it firmly so as to prevent the microscope from turning.

[G2] Star fixation knob

Lock the first arm by tightening it firmly so as to prevent the microscope from turning.

[G3] 5-pin plug/socket (manual focus without it)

Connect the foot switch.

[G4] Star fixation knob

Tighten it in order to lock the light source box to make it not revolve when moving or

storing.

[G5] Inner hexagonal tightening bolt

Lock the stand pillar and the light source box

[G6] 110/220V selecting switch.

[G7] Power plug

Input power supply

[G8] Star fixation knob (with sterilized cap)

Tighten it in order to lock the little arm and make the microscope not move vertically.




Fig.10

[E1] lamp-house box

[E3] table

[F2] upside of the column




Fig.11

[E2] clamp rack

[F1] spanner hole

[F3] underside of the column

The camera package includes the following components:

1. 1/4" CCD	1pc
2. CCD adapter	1pc
3. Beam splitter	1pc
4. transformer	1pc
5. Video frequency line (with adapter)	1pc
6. 2mm inner hexagram screw drive	1pc

★ Design or specifications are subject to change without prior notice.





Wellkang Ltd t/a Wellkang Tech Consulting,
Suite B,29 Harley Street
LONDON W1G 9QR
England, United Kingdom

MEMORANDUM



66 Vision Tech Co., Ltd.

Add: 34 Daru lane, Suzhou,
Jiangsu province P.R.C
Tel: (86) 512-67276129/67271504
Fax: (86) 512-67771975/65244789
Cable: 7729
Zip Code: 215005
Web Site: www.66vision.com

PACKING LIST

Content	Quan.
Stand pillar	1
Fixing bolt of column	1
Exponent	1
Base	1
Power supply wires	1
Foot control switch	1
Manipulating handle	2
Lens	1
Arm	1
Fiber optic	1
Sterilized cap of star fixation knob	4
Sterilized cap of manipulating handle	2
Sterilized cap of magnification knob	2
Medical halogen bulb with cool reflection (12v100w)	1
Fuse, T type AC250 1.25A、 T type AC125 2.5A	4
6mm inner hexagon spanner	1
10mm inner hexagon spanner	1
P54M bulb-socket	2

Spare parts:




Content	Quan.	Remark
1/4" CCD	1	
CCD adapter	1	
Beam splitter	1	
Transformer	1	
Video wires (with adapter)	1	
2mm inner hexagon spanner	1	

(With the pack sketch map on the back)

Suzhou 66 Vision Tech Co.,Ltd

YZ20P5 Operation Microscope instruction Manual

Pack Sketch Map



Sterilized cap , Fuse

Inner hexagon spanner , bulb-socket